


Anti-Cd96/Tactile Rabbit Monoclonal Antibody For Fc, Elisa - Q3U0X8

Item Number: CM0004167



Introduction

Rabbit monoclonal antibody against mouse CD96/TACTILE. Validated for flow cytometry (FC) and ELISA. Targets T-cell surface protein involved in NK cell adhesion. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD96/TACTILE Rabbit mAb
Remarks/Alias	Tactile; 1700109I12Rik
Species	Mouse
GeneID (Human)	10225
GeneID (Mouse)	84544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-536 of mouse CD96/TACTILE (NP_115854.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Mouse
SWISS	Q3U0X8
Protein Weight	67kDa
Shipping	Ice bag