

Anti-Cd97 (Adgre5) Rabbit Polyclonal Antibody For Wb And Elisa - P48960

Item Number: CM0028383



Introduction

High-quality rabbit polyclonal antibody targeting human CD97 (ADGRE5). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying leukocyte adhesion, migration, and inflammatory diseases.

[Learn More](#)

Parameter	Value
Product Name	CD97 Rabbit pAb
Remarks/Alias	CD97; TM7LN1
Species	Human
GeneID (Human)	976
GeneID	976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-435 of human CD97 (NP_001775.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48960
Protein Weight	92kDa
Shipping	Ice bag