

Anti-Cd97 () Monoclonal Antibody For Flow Cytometry - P48960

Item Number: CM0008164



Introduction

-conjugated anti-human CD97 monoclonal antibody validated for flow cytometry. Reacts with human CD97, a leukocyte adhesion GPCR involved in migration and inflammation. Suitable for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD97 mAb
Remarks/Alias	CD97; TM7LN1
Species	Human
Gene ID	976
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-398 of human CD97 (NP_510966.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48960
Protein Weight	92kDa
Shipping	Ice bag