


Anti-Cd97 Rabbit Monoclonal Antibody For Wb, Elisa - P48960

Item Number: CM0007119



Introduction

Anti-CD97 rabbit monoclonal antibody for WB and ELISA applications. Specific to human CD97, an adhesion GPCR involved in leukocyte migration and inflammatory responses. Predicted molecular weight: 92 kDa.

[Learn More](#)

Parameter	Value
Product Name	CD97 Rabbit mAb
Remarks/Alias	CD97; TM7LN1
Species	Human
GeneID (Human)	976
GeneID	976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-398 of human CD97 (NP_510966.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P48960
Protein Weight	92 kDa
Shipping	Ice bag