

Anti-Cd97 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Z0M6

Item Number: CM0020489

Introduction

Rabbit polyclonal antibody targeting mouse CD97 (Adgre5), validated for WB and ELISA. Suitable for studying leukocyte adhesion, migration, and GPCR signaling pathways in immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD97 Rabbit pAb
Remarks/Alias	Cd97; TM7LN1; CD97
Species	Mouse
Gene ID (Human)	976
Gene ID	26364
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of mouse CD97. (NP_036055.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Z0M6
Protein Weight	90kDa
Shipping	Ice bag