

Anti-Cd98 Monoclonal Antibody For Fc - P10852

Item Number: CM0009616



Introduction

High-quality rabbit monoclonal antibody against mouse CD98 (SLC3A2), validated for flow cytometry. Suitable for studying amino acid transport, integrin signaling, and tumor progression in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD98 mAb
Remarks/Alias	4F2; Cd98; Ly10; Mdu1; 4F2HC; Ly-10; NACAE; Ly-m10; Mgp-2hc
Species	Mouse
GeneID (Human)	6520
GeneID	17254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-526 of mouse CD98 (NP_032603.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P10852
Protein Weight	58kDa
Shipping	Ice bag