

Anti-Cd98 Rabbit Polyclonal Antibody For Wb - P10852

Item Number: CM0022170



Introduction

High-quality rabbit polyclonal antibody against mouse CD98 (SLC3A2) for Western blotting and ELISA. Validated in mouse and rat tissues. Recombinant immunogen, molecular weight 58 kDa. Ideal for amino acid transport and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD98 Rabbit pAb
Remarks/Alias	4F2; Cd98; Ly10; Mdu1; 4F2HC; Ly-10; NACAE; Ly-m10; Mgp-2hc; CD98
Species	Mouse
GeneID	17254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-526 of mouse CD98. (NP_032603.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10852
Protein Weight	58kDa
Shipping	Ice bag