

Anti-Cd98/Slc3A2 Monoclonal Antibody For Wb, If, Ihc, Elisa - P08195

Item Number: CM0007849



Introduction

Rabbit monoclonal antibody against human CD98/SLC3A2 (4F2hc), validated for WB, IF, IHC, ELISA. Reacts with human, mouse, rat. Suitable for amino acid transporter studies and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD98/SLC3A2 Rabbit mAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human
GeneID (Human)	6520
GeneID	6520
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of Human SLC3A2/CD98hc (NP_002385.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08195
Protein Weight	58-71 kDa
Shipping	Ice bag