


Anti-Cd98/Slc3A2 Monoclonal Antibody For Wb, If/Icc, Elisa - P08195 / P10852

Item Number: CM0002884



Introduction

Rabbit monoclonal antibody against human/mouse CD98/SLC3A2, validated for WB, IF/ICC, ELISA. Reacts with human and rat. Target is a heavy chain of amino acid transporters involved in amino acid exchange and integrin signaling.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated]CD98/SLC3A2 Rabbit PolymAb®
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human/Mouse
Gene ID (Human)	6520
Gene ID	6520/17254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of human/mouse SLC3A2/CD98hc (NP_002385.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P08195/P10852
Protein Weight	58kDa
Shipping	Ice bag