

Anti-Cd98/Slc3A2 Rabbit Monoclonal Antibody For Western Blot, Immunofluorescence, Elisa - P08195

Item Number: CM0006767



Introduction

High-quality rabbit monoclonal antibody against human CD98/SLC3A2 (58 kDa, P08195). Validated for WB, IF/ICC, and ELISA. Targets the heavy chain of amino acid transporters involved in nutrient uptake, integrin signaling, and viral entry. Ideal for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CD98/SLC3A2 Rabbit mAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human
GeneID (Human)	6520
GeneID	6520
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 461-630 of human SLC3A2/CD98hc (NP_002385.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P08195
Protein Weight	58kDa
Shipping	Ice bag