


Anti-Cd99 Polyclonal Antibody For Wb, Elisa - P14209

Item Number: CM0023366



Introduction

CD99 rabbit polyclonal antibody validated for western blot (WB) and ELISA. Targets human CD99 (MIC2, 19 kDa) involved in T-cell adhesion and leukocyte extravasation. Ideal for immunology and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CD99 Rabbit pAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X; CD99
Species	Human
GeneID (human)	4267
GeneID	4267
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 23-124 of human CD99 (NP_002405.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14209
Protein Weight	19kDa
Shipping	Ice bag