


Anti-Cd99 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14209

Item Number: CM0020880



Introduction

CD99 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human samples. Targets CD99 antigen involved in T-cell adhesion and leukocyte transmigration. Recombinant immunogen (aa 23-185).

[Learn More](#)

Parameter	Value
Product Name	CD99 Rabbit pAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X; CD99
Species	Human
Gene ID (Human)	4267
Gene ID	4267
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P14209
Protein Weight	19kDa
Shipping	Ice bag