


Anti-Cd99 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Fc, Elisa - P14209

Item Number: CM0010785



Introduction

High-quality CD99 rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, FC, ELISA. Reacts with human CD99 antigen (MIC2, HBA71) involved in T-cell adhesion. GeneID: 4267, SWISS: P14209.

[Learn More](#)

Parameter	Value
Product Name	CD99 Rabbit mAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X
Species	Human
GeneID (Human)	4267
GeneID	4267
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 24-123 of human CD99 (NP_002405.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P14209
Protein Weight	16-19kDa
Shipping	Ice bag