

Anti-CD99L2 Rabbit Polyclonal Antibody For Wb - Q8Tcz2

Item Number: CM0022037

Introduction

CamelBio offers anti-CD99L2 rabbit polyclonal antibody, validated for WB and ELISA. Reacts with human, mouse, and rat CD99L2. CD99L2 plays a role in leukocyte extravasation and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	CD99L2 Rabbit pAb
Remarks/Alias	CD99B; MIC2L1; CD99L2
Species	Human
Gene ID (Human)	83692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-185 of human CD99L2 (NP_113650.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCZ2
Protein Weight	28kDa
Shipping	Ice bag