

Anti-Cda Monoclonal Antibody For Wb, Ihc-P, Elisa - P32320

Item Number: CM0007509



Introduction

Rabbit monoclonal antibody against human CDA (cytidine deaminase), suitable for WB, IHC-P, ELISA. Validated in human and mouse species. Predicted molecular weight 16 kDa. Ideal for nucleotide metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	CDA Rabbit mAb
Remarks/Alias	CDD; CDA
Species	Human
GeneID (Human)	978
GeneID	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA(NP_001776.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32320
Protein Weight	16kDa
Shipping	Ice bag