

Anti-Cda Polyclonal Antibody For Wb, If/Icc, Elisa - P32320

Item Number: CM0024184



Introduction

Rabbit polyclonal antibody targeting cytidine deaminase (CDA), validated for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat CDA. Ideal for nucleoside metabolism studies. Accession: P32320.

[Learn More](#)

Parameter	Value
Product Name	CDA Rabbit pAb
Remarks/Alias	CDD, CDA
Species	Human
Gene ID	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA (NP_001776.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32320
Protein Weight	16 kDa
Shipping	Ice bag