

Anti-Cda Rabbit Monoclonal Antibody For Ihc-P, Fc (Intra), Elisa - P32320

Item Number: CM0005965



Introduction

CamelBio offers Anti-CDA Rabbit Monoclonal Antibody validated for IHC-P, FC (intra), ELISA. Reacts with human and mouse CDA, a cytidine deaminase involved in pyrimidine salvage. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	CDA Rabbit mAb
Remarks/Alias	CDD; CDA
Species	Human
GeneID (Human)	978
GeneID	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA (NP_001776.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse
SWISS	P32320
Protein Weight	16kDa
Shipping	Ice bag