

Anti-Cda Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P32320

Item Number: CM0015100



Introduction

Cytidine deaminase (CDA) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse CDA (16 kDa); ideal for nucleotide metabolism and granulocyte research.

[Learn More](#)

Parameter	Value
Product Name	CDA Rabbit pAb
Remarks/Alias	CDD; CDA
Species	Human
GeneID (Human)	978
GeneID	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA (NP_001776.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32320
Protein Weight	16kDa
Shipping	Ice bag