


Anti-Cdc123 Rabbit Polyclonal Antibody For Wb And Elisa - 075794

Item Number: CM0014549



Introduction

High-quality rabbit polyclonal antibody against human CDC123 (D123, C10orf7) validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for translation initiation studies and eIF2 complex research.

[Learn More](#)

Parameter	Value
Product Name	CDC123 Rabbit pAb
Remarks/Alias	D123; C10orf7; CDC123
Species	Human
GeneID	8872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human CDC123 (NP_006014.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075794
Protein Weight	39kDa
Shipping	Ice bag