


Anti-Cdc16 Rabbit Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q13042

Item Number: CM0012732



Introduction

High-quality anti-CDC16 rabbit monoclonal antibody validated for IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat CDC16. Suitable for cell cycle and ubiquitin ligase research.

[Learn More](#)

Parameter	Value
Product Name	CDC16 Rabbit mAb
Remarks/Alias	APC6; CUT9; ANAPC6; CDC16Hs; CDC16
Species	Human
Gene ID (Human)	8881
Gene ID	8881
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CDC16 (Q13042).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13042
Protein Weight	72kDa
Shipping	Ice bag