


Anti-Cdc16 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13042

Item Number: CM0025708



Introduction

High-quality rabbit polyclonal antibody targeting human CDC16 (APC6). Validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for studying APC/C complex function, cell cycle regulation, and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	CDC16 Rabbit pAb
Remarks/Alias	APC6; CUT9; ANAPC6; CDC16Hs; CDC16
Species	Human
GeneID (Human)	8881
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 371-620 of human CDC16 (NP_003894.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13042
Protein Weight	72kDa
Shipping	Ice bag