


Anti-Cdc20 Polyclonal Antibody For Wb, If/Icc, Elisa - Q12834

Item Number: CM0014786



Introduction

Rabbit polyclonal antibody against human CDC20, validated for WB, IF/ICC, and ELISA. Recognizes human and mouse CDC20, a critical APC/C activator and mitotic checkpoint component. Suitable for cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CDC20 Rabbit pAb
Remarks/Alias	CDC20A; OOMD14; p55CDC; OZEMA14; bA276H19.3; CDC20
Species	Human
Gene ID (Human)	991
Gene ID	991
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CDC20 (NP_001246.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12834
Protein Weight	55kDa
Shipping	Ice bag