


Anti-Cdc20 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - Q12834

Item Number: CM0015382



Introduction

High-quality rabbit polyclonal antibody targeting human CDC20 (p55CDC), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, rat. Ideal for cell cycle, mitosis, and APC/C research.

[Learn More](#)

Parameter	Value
Product Name	CDC20 Rabbit pAb
Remarks/Alias	CDC20A; OOMD14; p55CDC; OZEMA14; bA276H19.3; CDC20
Species	Human
GeneID	991
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CDC20 (NP_001246.2?).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12834
Protein Weight	55kDa
Shipping	Ice bag