


Anti-Cdc23 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ujx2

Item Number: CM0012589



Introduction

Recombinant rabbit monoclonal antibody targeting human CDC23 (APC8, ANAPC8), validated for WB and ELISA. Suitable for cell cycle and ubiquitin ligase research, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CDC23 Rabbit mAb
Remarks/Alias	APC8; CUT23; ANAPC8; CDC23
Species	Human
GeneID (Human)	8697
GeneID	8697
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 498-597 of human CDC23 (Q9UJX2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJX2
Protein Weight	69 kDa
Shipping	Ice bag