


Anti-Cdc23 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ujx2

Item Number: CM0017187



Introduction

Rabbit polyclonal antibody against human CDC23 (APC8). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for cell cycle and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	CDC23 Rabbit pAb
Remarks/Alias	APC8; CUT23; ANAPC8; CDC23
Species	Human
GeneID (Human)	8697
GeneID	8697
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human CDC23 (NP_004652.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJX2
Protein Weight	69kDa
Shipping	Ice bag