


Anti-Cdc25A Rabbit Monoclonal Antibody For Wb, Elisa - P30304

Item Number: CM0004684



Introduction

Rabbit monoclonal antibody targeting human CDC25A (M-phase inducer phosphatase 1). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CDC25A Rabbit mAb
Remarks/Alias	CDC25A2
Species	Human
GeneID (Human)	993
GeneID	993
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 376-524 of human CDC25A (NP_001780.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30304
Protein Weight	54kDa/59kDa
Shipping	Ice bag