


Anti-Cdc25A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30304

Item Number: CM0017697



Introduction

Rabbit anti-CDC25A polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat CDC25A (59 kDa). Key cell cycle phosphatase for mitosis research.

[Learn More](#)

Parameter	Value
Product Name	CDC25A Rabbit pAb
Remarks/Alias	CDC25A2; CDC25A
Species	Human
GeneID (Human)	993
GeneID	993
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30304
Protein Weight	59kDa
Shipping	Ice bag