

Anti-Cdc25B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P30305

Item Number: CM0005436

Introduction

Rabbit monoclonal antibody against human CDC25B, a dual specificity phosphatase involved in cell cycle regulation. Validated for WB, IHC-P, ELISA. Ideal for studies on mitotic progression and cell division.

[Learn More](#)

Parameter	Value
Product Name	CDC25B Rabbit mAb
Remarks/Alias	CDC25B; cell division cycle 25B
Species	Human
GeneID (Human)	994
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDC25B (P30305).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P30305
Protein Weight	65kDa
Shipping	Ice bag