

## Anti-Cdc25C Monoclonal Antibody For Wb, If/Icc, Elisa - P30307

Item Number: CM0002137



### Introduction

Rabbit monoclonal anti-Cdc25C antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat Cdc25C, a critical mitotic regulator and CDK1 phosphatase. Ideal for cell cycle and oncology research.

[Learn More](#)

| Specification    | Details                                                                                            |
|------------------|----------------------------------------------------------------------------------------------------|
| Product Name     | Cdc25C Rabbit mAb                                                                                  |
| Remarks/Alias    | CDC25; PPP1R60; Cdc25C                                                                             |
| Species          | Human                                                                                              |
| GeneID (Human)   | 995                                                                                                |
| GeneID           | 995                                                                                                |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cdc25C (P30307). |
| Source           | Rabbit                                                                                             |
| Category         | Monoclonal Antibodies                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                  |
| SWISS            | P30307                                                                                             |
| Protein Weight   | 53kDa                                                                                              |
| Shipping         | Ice bag                                                                                            |