

Anti-Cdc25C Rabbit Polyclonal Antibody For Wb, Elisa - P30307

Item Number: CM0022399



Introduction

High-quality CDC25C rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat CDC25C. Ideal for cell cycle and mitosis research. SWISS: P30307.

[Learn More](#)

Parameter	Value
Product Name	CDC25C Rabbit pAb
Remarks/Alias	CDC25; PPP1R60; CDC25C
Species	Human
GeneID (Human)	995
GeneID	995
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CDC25C (NP_001781.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30307
Protein Weight	53kDa
Shipping	Ice bag