

Anti-Cdc34 Polyclonal Antibody For Wb, If/Icc, Elisa - P49427

Item Number: CM0028900



Introduction

Rabbit polyclonal antibody against human CDC34 (UBE2R1, E2-CDC34). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for ubiquitination and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CDC34 Rabbit pAb
Remarks/Alias	UBC3; UBCH3; UBE2R1; E2-CDC34; CDC34
Species	Human
Gene ID	997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-236 of human CDC34 (NP_004350.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49427
Protein Weight	27kDa
Shipping	Ice bag