

Anti-Cdc34 Rabbit Monoclonal Antibody For Wb, Elisa - P49427

Item Number: CM0012630



Introduction

Rabbit monoclonal antibody against human CDC34 (UBE2R1), validated for WB and ELISA. Recognizes the ubiquitin-conjugating enzyme E2 R1, essential for ubiquitin-mediated proteolysis and cell cycle regulation. Ideal for research on protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	CDC34 Rabbit mAb
Remarks/Alias	UBC3; UBCH3; UBE2R1; E2-CDC34; CDC34
Species	Human
GeneID (Human)	997
GeneID	997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-236 of human CDC34 (P49427).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49427
Protein Weight	27 kDa
Shipping	Ice bag