


Anti-Cdc37 Rabbit Monoclonal Antibody For Wb And Elisa - Q16543

Item Number: CM0010140



Introduction

High-quality rabbit monoclonal antibody targeting CDC37 (Hsp90 co-chaperone). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Useful for kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	CDC37 Rabbit mAb
Remarks/Alias	P50CDC37; CDC37
Species	Human
GeneID (Human)	11140
GeneID	11140
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDC37 (Q16543)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16543
Protein Weight	44kDa
Shipping	Ice bag