


Anti-Cdc37 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q16543

Item Number: CM0026964



Introduction

Anti-CDC37 Rabbit Polyclonal Antibody validated for WB, IP, ELISA in Human and Mouse samples. Targets Hsp90 co-chaperone Cdc37 (44kDa), involved in kinase stabilization. Ideal for chaperone biology research.

[Learn More](#)

Parameter	Value
Product Name	CDC37 Rabbit pAb
Remarks/Alias	P50CDC37; CDC37
Species	Human
GeneID (Human)	11140
GeneID	11140
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human CDC37 (NP_008996.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16543
Protein Weight	44kDa
Shipping	Ice bag