

Anti-Cdc40 Rabbit Monoclonal Antibody (Prp17/Prpf17) For Wb, Ihc-P, Elisa - O60508

Item Number: CM0002617



Introduction

High-quality rabbit monoclonal antibody targeting human CDC40 (PRP17, EHB3). Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for spliceosome research and embryonic brain development studies.

[Learn More](#)

Parameter	Value
Product Name	CDC40 Rabbit mAb
Remarks/Alias	EHB3; PCH15; PRP17; PRPF17; CDC40
Species	Human
Gene ID (Human)	51362
Gene ID	51362
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60508
Protein Weight	66 kDa
Shipping	Ice bag