

Anti-Cdc42 Rabbit Monoclonal Antibody For Wb And Elisa - P60953

Item Number: CM0003769



Introduction

This rabbit monoclonal antibody targets human CDC42 (also mouse, rat) and is validated for Western blot and ELISA. Ideal for studying cell division control, actin cytoskeleton, and epithelial polarization.

[Learn More](#)

Parameter	Value
Product Name	CDC42 Rabbit mAb
Remarks/Alias	TKS; G25K; CDC42Hs; CDC42
Species	Human
GeneID (Human)	998
GeneID	998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-191 of human CDC42 (P60953)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60953
Protein Weight	21kDa
Shipping	Ice bag