

Anti-Cdc42 Rabbit Polyclonal Antibody For If/Icc, Elisa - P60953

Item Number: CM0016895



Introduction

Rabbit polyclonal antibody against human CDC42, optimized for IF/ICC and ELISA. Targets a small GTPase involved in cell polarization, migration, and cytoskeletal remodeling. Suitable for cellular biology research.

[Learn More](#)

Parameter	Value
Product Name	CDC42 Rabbit pAb
Remarks/Alias	TKS; G25K; CDC42Hs; CDC42
Species	Human
GeneID (human)	998
GeneID	998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-186 of human CDC42 (NP_001782.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P60953
Protein Weight	21kDa
Shipping	Ice bag