

Anti-Cdc42Ep1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q00587

Item Number: CM0002799

Introduction

Rabbit monoclonal antibody against human CDC42EP1 (BORG5/MSE55), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for actin cytoskeleton and endothelial cell research.

[Learn More](#)

Parameter	Value
Product Name	CDC42EP1 Rabbit mAb
Remarks/Alias	CEP1; BORG5; MSE55; CDC42EP1
Species	Human
GenelD (human)	11135
GenelD	11135
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 292-391 of human CDC42EP1 (Q00587).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00587
Protein Weight	40kDa
Shipping	Ice bag