

Anti-Cdc45 Polyclonal Antibody For Wb, Elisa - 075419

Item Number: CM0022377



Introduction

Rabbit polyclonal antibody targeting human CDC45, validated for WB and ELISA. Reacts with human, mouse, rat. Essential core component of the CMG helicase for DNA replication initiation. Ideal for cell cycle and chromosome research.

[Learn More](#)

Parameter	Value
Product Name	CDC45 Rabbit pAb
Remarks/Alias	CDC45L; MGORS7; CDC45L2; PORC-PI-1; CDC45
Species	Human
Gene ID (Human)	8318
Gene ID	8318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CDC45 (NP_003495.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75419
Protein Weight	66kDa
Shipping	Ice bag