

Anti-Cdc5L Rabbit Polyclonal Antibody For Wb, Elisa - Q99459

Item Number: CM0025828



Introduction

Rabbit polyclonal antibody targeting human CDC5L (aa 563-802) for WB and ELISA applications. Cross-reacts with human and mouse. Ideal for cell cycle and pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	CDC5L Rabbit pAb
Remarks/Alias	CDC5; CEF1; PCDC5RP; CDC5-LIKE; dj319D22.1; CDC5L
Species	Human
Gene ID (Human)	988
Gene ID	988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 563-802 of human CDC5L (NP_001244.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99459
Protein Weight	92kDa
Shipping	Ice bag