


Anti-Cdc6 Polyclonal Antibody For Wb - Q99741

Item Number: CM0021434



Introduction

Rabbit polyclonal antibody against human Cdc6, validated for Western blot and ELISA. Useful for DNA replication and cell cycle checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	Cdc6 Rabbit pAb
Remarks/Alias	CDC18L; HsCDC6; MGORS5; HsCDC18; Cdc6
Species	Human
GeneID (Human)	990
GeneID	990
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-560 of human Cdc6 (Q99741)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99741
Protein Weight	63kDa
Shipping	Ice bag