

Anti-Cdca3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99618

Item Number: CM0004464

Introduction

High-quality anti-CDCA3 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human CDCA3, a cell division cycle-associated protein involved in mitotic entry. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CDCA3 Rabbit mAb
Remarks/Alias	GRCC8; TOME1; TOME-1; CDCA3
Species	Human
Gene ID	83461 (Human)
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDCA3 (Q99618).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q99618
Protein Weight	29kDa
Shipping	Ice bag