

Anti-Cdca3 Rabbit Polyclonal Antibody For Wb, Elisa - Q99618

Item Number: CM0023838

Introduction

Anti-CDCA3 polyclonal antibody from rabbit, validated for Western blot (WB) and ELISA. Detects human and mouse CDCA3, a mitotic entry trigger. Ideal for cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CDCA3 Rabbit pAb
Remarks/Alias	GRCC8; TOME1; TOME-1; CDCA3
Species	Human
GeneID (Human)	83461
GeneID	83461
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human CDCA3 (NP_001284531.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99618
Protein Weight	29kDa
Shipping	Ice bag