

Anti-Cdca5 Polyclonal Antibody For Wb - Q96Ff9

Item Number: CM0022910

Introduction

CDCA5 Rabbit pAb is a polyclonal antibody raised in rabbit against human CDCA5 (sororin). Validated for WB and ELISA, it reacts with human, mouse, and rat samples. Ideal for studying sister chromatid cohesion and cell division.

[Learn More](#)

Parameter	Value
Product Name	CDCA5 Rabbit pAb
Remarks/Alias	SORORIN; CDCA5
Species	Human
Gene ID	113130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-150 of human CDCA5 (NP_542399.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FF9
Protein Weight	28kDa
Shipping	Ice bag