


Anti-Cdca8 Rabbit Polyclonal Antibody For Wb, Elisa - Q53HL2

Item Number: CM0015868



Introduction

Rabbit polyclonal antibody against human CDCA8, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Validated for research on mitosis and cell division.

[Learn More](#)

Parameter	Value
Product Name	CDCA8 Rabbit pAb
Remarks/Alias	BOR; DasraB; MESRGP; BOREALIN; CDCA8
Species	Human
Gene ID (Human)	55143
Gene ID	55143
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-124 of human CDCA8 (NP_060571.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53HL2
Protein Weight	31kDa
Shipping	Ice bag