


## Anti-Cdca8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q53HI2

Item Number: CM0013098



### Introduction

Rabbit polyclonal antibody targeting CDCA8 (Borealin), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mitotic regulation and cell cycle research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CDCA8 Rabbit pAb                                                                                                 |
| Remarks/Alias    | BOR; DasraB; MESRGP; BOREALIN; CDCA8                                                                             |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 55143                                                                                                            |
| GeneID           | 55143                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CDCA8 (NP_060571.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q53HL2                                                                                                           |
| Protein Weight   | 31kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |