

Anti-Cdcp1 Monoclonal Antibody For Fc - Q9H5V8

Item Number: CM0010225

Introduction

-conjugated rabbit monoclonal antibody targeting human CDCP1 (CD318, TRASK). For flow cytometry. CDCP1 is a cell adhesion protein involved in tumor progression and metastasis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CDCP1 mAb
Remarks/Alias	CD318; TRASK; SIMA135
Species	Human
GeneID (Human)	64866
GeneID	64866
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-666 of human CDCP1 (NP_073753.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9H5V8
Protein Weight	37kDa/72kDa/92kDa
Shipping	Ice bag