

Anti-Cdcp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H5V8

Item Number: CM0028278

Introduction

CDCP1 rabbit polyclonal antibody for WB and ELISA. Reacts with human and mouse CDCP1, a cell adhesion and tumor progression marker with roles in leukemia diagnosis and stem cell subsets.

[Learn More](#)

Parameter	Value
Product Name	CDCP1 Rabbit pAb
Remarks/Alias	CD318; TRASK; SIMA135; CDCP1
Species	Human
GeneID (Human)	64866
GeneID	64866
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-343 of human CDCP1 (NP_835488.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H5V8
Protein Weight	93kDa
Shipping	Ice bag