

Anti-Cdh11 Polyclonal Antibody For Wb, Elisa - P55287

Item Number: CM0029980

Introduction

Explore cell adhesion and migration with this rabbit anti-CDH11 polyclonal antibody. Validated for WB and ELISA, it recognizes human CDH11 (OB-cadherin) and cross-reacts with mouse and rat. Suitable for cadherin and focal adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CDH11 Rabbit pAb
Remarks/Alias	OB, ESWS, CAD11, CDHOB, OSF-4, TBHS2, CDH11
Species	Human
Gene ID (Human)	1009
Gene ID	1009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-400 of human CDH11 (NP_001788.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55287
Protein Weight	88kDa
Shipping	Ice bag