

## Anti-Cdh11 Polyclonal Antibody For Wb, If/Icc, Elisa - P55287

Item Number: CM0013118

### Introduction

Anti-CDH11 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human and mouse. Targets cadherin-11, a cell adhesion protein involved in focal adhesion assembly and migration, relevant for bone and brain research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CDH11 Rabbit pAb                                                                                                   |
| Remarks/Alias    | OB; ESWS; CAD11; CDHOB; OSF-4; TBHS2; CDH11                                                                        |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 1009                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-400 of human CDH11 (NP_001788.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | P55287                                                                                                             |
| Protein Weight   | 88kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |